

19



Europäisches Patentamt
European Patent Office
Office européen des brevets

11

Publication number:

0 182 546
A2

D5

12

EUROPEAN PATENT APPLICATION

21

Application number: 85308070.3

51

Int. Cl.: **G 01 N 33/541**

22

Date of filing: 06.11.85

30

Priority: 12.11.84 GB 8428487

71

Applicant: **POLYCLONAL ANTIBODIES LIMITED,**
14/15 Newbury Street, London EC1A 7HU (GB)

43

Date of publication of application: 28.05.86
Bulletin 86/22

72

Inventor: Landon, John, 85 West Street,
Harrow-on-the-Hill Middlesex (GB)

84

Designated Contracting States: **BE CH DE FR GB IT LI**
NL SE

74

Representative: **Wain, Christopher Paul et al, A.A.**
THORNTON & CO. Northumberland House 303-306 High
Holborn, London WC1V 7LE (GB)

54

Immunoassays.

57

In the immunoassay of an antigen using the second antibody technique, the second antibody is raised against the Fc component only of the first antibody. This reduces, or reduces the risk of, interference in the primary reaction by the second antibody. In a modified «two-site» immunoassay, a second antibody is used, directed against the Fc component of the unlabelled first antibody, to insolubilise the first antibody. Single reagent immunoassays may be effected using the modified second antibody technique of the invention.

EP 0 182 546 A2

R₅

⑫

EUROPEAN PATENT APPLICATION

⑰ Application number: 85308070.3

⑤① Int. Cl.⁴: **G 01 N 33/541**

⑰ Date of filing: 06.11.85

③① Priority: 12.11.84 GB 8428487

⑦① Applicant: **POLYCLONAL ANTIBODIES LIMITED,**
14/15 Newbury Street, London EC1A 7HU (GB)

④③ Date of publication of application: 28.05.86
Bulletin 86/22

⑦② Inventor: **Landon, John, 85 West Street,**
Harrow-on-the-Hill Middlesex (GB)

⑧④ Designated Contracting States: **BE CH DE FR GB IT LI**
NL SE

⑧⑧ Date of deferred publication of search
report: 09.09.87 Bulletin 87/37

⑦④ Representative: **Wain, Christopher Paul et al, A.A.**
THORNTON & CO. Northumberland House 303-306 High
Holborn, London WC1V 7LE (GB)

⑤④ **Immunoassays.**

⑤⑦ In the immunoassay of an antigen using the second antibody technique, the second antibody is raised against the Fc component only of the first antibody. This reduces, or reduces the risk of, interference in the primary reaction by the second antibody. In a modified "two-site" immunoassay, a second antibody is used, directed against the Fc component of the unlabelled first antibody, to insolubilise the first antibody. Single reagent immunoassays may be effected using the modified second antibody technique of the invention.

EP 0 182 546 A3



European Patent
Office

EUROPEAN SEARCH REPORT

0182546

Application number

EP 85 30 8070


DOCUMENTS CONSIDERED TO BE RELEVANT			
Category	Citation of document with indication, where appropriate, of relevant passages	Relevant to claim	CLASSIFICATION OF THE APPLICATION (Int. Cl. 4)
X	US-A-4 481 298 (R.O. CONE et al.) * complete *	1-10	G 01 N 33/541 G 01 N 33/563// G 01 N 33/74
A	US-A-4 298 592 (W.H.T. LIN et al.) * column 1, line 57 - column 2, line 57 *	1,2,5	
A	US-A-4 469 787 (J.W. WOODS et al.) * abstract *	1,2,5	
A,P	US-A-4 548 908 (N. KAMEDA) * abstract *	1	
A,P	US-A-4 522 922 (J. CARRO et al.) * claims *	1	TECHNICAL FIELDS SEARCHED (Int. Cl. 4) G 01 N 33/00
The present search report has been drawn up for all claims			
Place of search BERLIN		Date of completion of the search 15-06-1987	Examiner GREEN C.H.
CATEGORY OF CITED DOCUMENTS			
X : particularly relevant if taken alone Y : particularly relevant if combined with another document of the same category A : technical background O : non-written disclosure P : intermediate document		T : theory or principle underlying the invention E : earlier patent document, but published on, or after the filing date D : document cited in the application L : document cited for other reasons & : member of the same patent family, corresponding document	

- 1 -

IMMUNOASSAYS

This invention relates to immunoassays.

Immunoassays depend upon the reversible and non-covalent binding of an antigen (Ag), which is the analyte to be measured by its specific antibodies (Ab), to form a complex (Ag:Ab), in accordance with the reaction:



At equilibrium, some of the free reactants will be combining to form more complex, and some of the complex will be dissociating to give the free reactants.

In order to enable, or improve the accuracy of, quantitation in immunoassays, it is known to use a labelled reagent. Among the labels commonly used are radioisotopes, enzymes and fluorophors.


In many immunoassays, it is necessary to separate the bound fraction from the free fraction once equilibrium has been established. A large number of physicochemical, immunological and other means of separating the antibody-bound and free fractions have been developed. Among these is the so-called second antibody approach. In accordance with this technique, there is included in the reaction mixture, in addition to the antibody (first antibody) directed against the antigen, another antibody (second

In one aspect, the present invention thus provides a method of immunoassay of an antigen wherein the antigen is reacted with a first antibody to form an antigen:antibody complex, and wherein a second antibody
 5 is used to bind to the first antibody to form a precipitate therewith, the precipitate being separated out, characterised in that the second antibody is directed only against the Fc component of the first antibody.

The method of the invention is applicable to all
 10 second antibody assays, including those using a labelled reagent, in which there is a separation step.


One particular form of immunoassay with which the present invention may be utilised is that known as a "two-site" technique for the assay of large molecules with
 15 more than one antigenic determinant. In a conventional "two-site" assay, the sample under assay is first incubated with an excess of an antibody covalently linked to an insoluble support and directed against one antigenic determinant (SP-Ab₁)

20



After careful washing, excess of a second, isotopically
 25 labelled antibody (Ab₂^{*}) directed against another part of the antigen molecule is added

30



Following further washes, the counts (in the case of a radiolabel) present in the solid phase fraction are determined and are related directly to the concentration of antigen present in the samples and standards.

the assay result. The so-called "single reagent" is, of course, a mixture of all the various reaction components required for the assay (excepting the sample or standard). In theory, single reagent label immuno-

5 assays should significantly reduce the complexity of label assays, and increase the speed and precision thereof. Also, they should reduce reagent manufacturing costs and also improve the stability of some isotopically labelled antigens.

10 However, despite these potential advantages, there have been only a very few proposals hitherto for single reagent label immunoassays and, so far as we are aware, none of these has been very satisfactory. For example, one commercially available kit for the radio-

15 immunoassay of thyroxine (T_4) used a single reagent comprising ^{125}I - T_4 , mixed with anti- T_4 antiserum and second antibody raised against the entire immunoglobulin G (IgG) of the species used to make the antiserum. However, in conducting the assay, the mixture of sample (or standard)

20 and single reagent had to be incubated for 30 minutes at $56^\circ C$, which was not very convenient. Furthermore, the assay gave relatively poor precision (up to 10% coefficients of variation) and high recovery (up to 125%) as compared with a conventional double antibody assay for T_4 .

25 In another proposed radioimmunoassay (RIA) for T_4 , it was found that the assay was accurate only if the "single reagent" was made up immediately before use. Again, this is not very convenient and somewhat negates the advantages theoretically achievable by "single reagent" procedures.

30 We have now found, in accordance with a further aspect of the present invention, by using a second antibody directed against the Fc fragment only of the first antibody, "single reagents" can be made up which are stable and reliable in subsequent use in an assay. Thus,

35 the invention includes a single reagent label immunoassay

with fluorescein isothiocyanate according to standard methods. Antiserum to HSA was raised in a sheep by standard procedures. Antisera to sheep IgG and to sheep Fc fragment were raised in rabbits by standard procedures.

- 5 Assay diluent buffer was sodium phosphate (50 mmol/l pH 7.4) containing bovine serum albumin (5 g/l) and sodium azide (1 g/l). Standard solutions of HSA were prepared in diluent buffer.

- 10 Experiments were performed using polystyrene test tubes which also served as cuvettes for measurement of fluorescence in a Perkin-Elmer Model 1000M filter fluorimeter (Kamel et al. (1980) Clin. Chem. 26, 1281-1284).

- For the fluoroimmunoassays with conventional reagent addition sequence, working tracer reagent was
15 prepared by diluting the fluorescein-labelled HSA to 2.6 mg/l (based on its protein content) in diluent buffer and working antiserum reagent by diluting the anti-HSA serum by 1,430-fold in diluent buffer.

- The assay procedure with conventional second-
20 antibody separation was as follows. To 50 µl of standard was added 200 µl of working tracer reagent followed by 100 µl of working antiserum reagent. After incubation for 1 h at room temperature, 50 µl of normal sheep serum (diluted 1,000-fold in diluent buffer) and 50 µl of second antiserum
25 (either anti-IgG or anti-Fc, each diluted 8-fold in diluent buffer) was added. After further incubation for 2 h at room temperature, the assay mixtures were centrifuged, the supernatants aspirated to waste, and the precipitates re-dissolved by the addition of 1.5 ml of a mixture of
30 methanol, water and 2 mol/l sodium hydroxide (50:50:1 parts by volume). The fluorescence of the re-dissolved precipitates (antibody-bound fraction of the labelled HSA) was measured by placing each test tube into the fluorimeter.

- 35 The fluoroimmunoassay standard curves obtained

2. Example of a single reagent immunoradiometric assay
employing anti-Fc second antibody: liquid-phase
two-site assay of human prolactin

Antisera to human prolactin were raised in
5 sheep and rabbits by standard procedures. Specific anti-
prolactin immunoglobulins were isolated from the sheep
antiserum by immunoaffinity purification and were
radioiodinated with ^{125}I by conventional procedures.
Antisera to rabbit IgG and to rabbit Fc fragment were
10 raised in sheep by standard procedures. Assay diluent
buffer was sodium phosphate (64 mmol/l, pH 7.4)
containing polyethylene glycol 6000 (30 g/l), bovine serum
albumin (5 g/l), normal rabbit serum (10 ml/l), normal
sheep serum (20 ml/l), and thiomersal (100 mg/l). Standard
15 solutions of human prolactin were prepared in horse serum.

The procedure for conventional liquid-phase
two-site assay was as follows. To 100 μl of standard
was added 100 μl of rabbit anti-prolactin serum (diluted
400-fold in diluent buffer) followed by 100 μl (approx.
20 100,000 counts/min) of ^{125}I -labelled sheep anti-prolactin
immunoglobulins. After incubation for 3 h at room
temperature, 100 μl of sheep anti-(rabbit (IgG) serum
(appropriately) diluted in 64 mmol/l sodium phosphate
buffer, pH 7.4) was added and after incubation for a
25 further 1 h at room temperature the assay mixtures were
centrifuged, the supernatants aspirated to waste and the
radioactivity in the precipitates counted. The standard
curve obtained is shown in Figure 3.

When the above experiment was repeated employing
30 an appropriate dilution of sheep anti-(rabbit Fc) second
antiserum in place of the anti-(rabbit IgG) serum, an
identical standard curve was obtained.

The experiment of Figure 3 was repeated employing
a mixture of the two specific anti-prolactin reagents that
35 had been prepared several days previously. To 100 μl of
standard was added 200 μl _____

0182546

- 11 -

conventional procedure - avoiding two pipetting steps
and the need for a second incubation, and with all
reagent additions completed at one time.

6. A single reagent label immunoassay for an antigen, in which the antigen is mixed with a single reagent mixture and, optionally after a separation step, the label is assayed, characterised in that the single
5 reagent mixture comprises a first antibody against the antigen and a second antibody against the Fc component of the first antibody.
7. An assay according to claim 6, wherein said
10 antigen has at least two antigenic determinants, and wherein said "single reagent" mixture comprises antibody raised in a first species against a determinant of the antigen under assay; labelled antibody raised in a second species against another determinant of the antigen;
15 and antibody raised in said second species against the Fc fragment of the IgG of the first species.
8. An immunoassay according to claim 6 or 7, wherein the label is a radioisotope, enzyme or fluorophore.
20
9. A "single reagent" mixture for use in the assay of claim 6, which comprises a first antibody against an antigen, and a second antibody against the Fc component of the first antibody.
25
10. A "single reagent" mixture for use in the assay of claim 7, which comprises antibody of a first animal species against an antigen, labelled antibody of a second animal species against the antigen, and antibody raised
30 in said first species against the Fc component of the IgG of the first species.

1/4




Fig. 1.

2/4




FIG. 2.

3/4




FIG. 3.

4/4




FIG. 4.